

# MiniPet

## FEATURES

- > Compact and engineered for precision
- > Accuracy and affordability for routine applications
- > Custom colors with your name and logo
- > Optional DeTipper Bottle
- > Uses standard 0—200  $\mu$ L Yellow Tips

## CUSTOM OPTIONS


- > Custom options available for OEM customers only
  - Dispense Volume: 5  $\mu$ L - 200  $\mu$ L
  - Unit Color
  - Cap Pad Printing
  - Packaging

## OEM REQUIREMENTS

- > A volume of 1000 units allows custom options.

## BASE MODEL

5 $\mu$ L	40 $\mu$ L
10 $\mu$ L	50 $\mu$ L
20 $\mu$ L	75 $\mu$ L
30 $\mu$ L	100 $\mu$ L
35 $\mu$ L	200 $\mu$ L



(Shown with optional tip)

## SPECIFICATIONS

Model	5 µL Cat# 987	10 µL Cat# 871	20 µL Cat# 872	25 µL Cat# 873	30 µL Cat# 916	35 µL Cat# 874	40 µL Cat# 875	50 µL Cat# 876	75 µL Cat# 926	100 µL Cat# 877	200 µL Cat# 918
Color	Red	White	Pink	Orange	Plum	Light Blue	Mint Green	Yellow	Green	Blue	Teal
Tips	Generic tips can be used with most MiniPets										
Liquid Contact Materials :	Polypropylene (disposable tip)										

## FUNCTIONAL PARAMETERS


Precision (full stroke)	2% CV or better		
Accuracy (full stroke)	±3%		
Performance (with generic Eppendorf style 0-200 µL Tip)	Volume	± % Accuracy	% Precision (CV)
	200 µL (full volume)	3	2
	100 µL (full volume)	3	2
	75 µL	3	2
	50 µL (full volume)	3	2
	35 µL- 45 µL	4	2
	20 µL- 30 µL	5	3
	10 µL, 11 µL and 14 µL	10	4
	5 µL	20	5
Dimensions (instrument only)			
Length	3.75 in. [9.53 cm]		
Handle diameter	.540 in. [1.37 cm]		
Weight	.2 oz [6.25 gm]		

## CUSTOM OPTIONS

Unit Color					
White <i>Omn# UN0001</i>	Red <i>PMS# 201C</i>	Dark Blue <i>PMS# 293C</i>	Yellow <i>PMS# 129C</i>	Mint Green <i>PMS# 346C</i>	Plum <i>PMS# 2593C</i>
Gray <i>PMS# 8C</i>	Pink <i>PMS# 230C</i>	Bright Yellow <i>PMS# 012C</i>	Green <i>PMS# 347C</i>	Red <i>PMS# 200C</i>	Aqua <i>PMS# 313U</i>
Tan <i>PMS# 451C</i>	Light Blue <i>PMS# 297C</i>	Medium Blue <i>PMS# 285C</i>	Purple <i>PMS# 2665C</i>	Orange <i>PMS# 1575C</i>	Teal <i>PMS# 3145C</i>

## PAD PRINTING OPTIONS

Pad Printing Ink Color	White or Black
------------------------	----------------



## MINIMUM OEM ORDER / SHIPMENT REQUIREMENTS

1,000 each

One time set up charge may apply.

PRECISION DISPENSING MICROPIPETTE FOR OEMS WORLDWIDE

[tricontinent.com](http://tricontinent.com)



**TriContinent Scientific, Inc.**  
12555 Loma Rica Drive  
Grass Valley, CA 95945  
USA  
Tel: 800-937-4738  
Fax: 530-273-2586  
[liquidhandling@tricontinent.com](mailto:liquidhandling@tricontinent.com)

**Gardner Denver Thomas GmbH**  
Livry-Gargan-Str. 10  
82256 Fürstenfeldbruck  
Germany  
Tel: +49 8141 2280 0  
Fax: +49 8141 8892136  
[thomas.de@gardnerdenver.com](mailto:thomas.de@gardnerdenver.com)

**Gardner Denver Hong Kong, Limited**  
Unit 1317-1318 Delta House, 3 On Yiu Street  
Siu Lek Yuen, Shatin, New Territories  
Hong Kong  
Tel: +852 26903502  
Fax: +852 27924598  
[thomas.hk@gardnerdenver.com](mailto:thomas.hk@gardnerdenver.com)

**TRICONTINENT**  
by Gardner Denver

The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as an aid to help in the selection of TriContinent products. It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability in connection there with. TriContinent does not warrant, guarantee or assume any obligation or liability in connection with this information.

Photos of products pictured in this catalog do not necessarily represent a specific model number. To obtain further information for custom options, contact your local TriContinent office.

Printed in USA Form No. MKT90006 J © Gardner Denver TriContinent. All rights reserved.